

Product Name: Cytokeratin 18 (6Y2) Rabbit Monoclonal Antibody
Catalog #: AMRe09735

Summary

Production Name	Cytokeratin 18 (6Y2) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	KRT18
Alternative Names	
Gene ID	3875.0
SwissProt ID	P05783.A synthetic peptide of human Cytokeratin 18

Application

Dilution Ratio	WB: 1:2000
Molecular Weight	48kDa


Background

Product Name: Cytokeratin 18 (6Y2) Rabbit Monoclonal Antibody
Catalog #: AMRe09735



Research Area

Image Data



Western blot detection of Cytokeratin 18 in K562, HeLa cell lysates using Cytokeratin 18 antibody (1:1000 diluted).

Note

For research use only.