

Product Name: D4 GDI (13R12) Rabbit Monoclonal Antibody
Catalog #: AMRe09769



Summary

Production Name	D4 GDI (13R12) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	ARHGDIB
Alternative Names	
Gene ID	397.0
SwissProt ID	P52566.A synthetic peptide of human D4 GDI

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	23kDa


Background

Product Name: D4 GDI (13R12) Rabbit Monoclonal Antibody
Catalog #: AMRe09769



Research Area

Image Data



Western blot detection of D4 GDI in Jurkat,U937 cell lysates using D4 GDI antibody(1:1000 diluted).

Note

For research use only.