

Product Name: DAZL (11G2) Rabbit Monoclonal Antibody
Catalog #: AMRe09800

Summary

Production Name	DAZL (11G2) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	DAZL
Alternative Names	DAZ homolog; DAZ like autosomal; DAZH; DAZL; DAZL1; DAZLA; Deleted in azoospermia like 1; Germline specific RNA binding protein; SPGY like autosomal; SPGYLA; Tpx2;
Gene ID	1618.0
SwissProt ID	Q92904.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	33kDa


Product Name: DAZL (11G2) Rabbit Monoclonal Antibody
Catalog #: AMRe09800

Background

RNA-binding protein, which is essential for gametogenesis. Plays a central role during spermatogenesis. May act by binding to the 3'-UTR of mRNA and thereby regulating the translation of key transcripts. RNA-binding protein, which is essential for gametogenesis in both males and females. Plays a central role during spermatogenesis. Acts by binding to the 3'-UTR of mRNA, specifically recognizing GUU triplets, and thereby regulating the translation of key transcripts (By similarity).

Research Area

Image Data



Western blot analysis of DAZL expression in human testis lysate.

Note

For research use only.