

Product Name: DDX21 (19Y4) Rabbit Monoclonal Antibody
Catalog #: AMRe09877

Summary

Production Name	DDX21 (19Y4) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	DDX21
Alternative Names	
Gene ID	9188.0
SwissProt ID	Q9NR30.A synthetic peptide of human DDX21

Application


Dilution Ratio	WB: 1:1000
Molecular Weight	87kDa

Background

Product Name: DDX21 (19Y4) Rabbit Monoclonal Antibody
Catalog #: AMRe09877

Research Area

Image Data



Western blot detection of DDX21 in HeLa cell lysates using DDX21 antibody(1:1000 diluted).

Note

For research use only.