

Product Name: DGKA (5X13) Rabbit Monoclonal Antibody
Catalog #: AMRe09943

Summary

Production Name	DGKA (5X13) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% New type preservative N and 0.05% BSA
Purification	Affinity purification

Immunogen

Gene Name	DGKA
Alternative Names	
Gene ID	1606.0
SwissProt ID	P23743.

Application


Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100
Molecular Weight	83 kDa

Product Name: DGKA (5X13) Rabbit Monoclonal Antibody
Catalog #: AMRe09943


Background

Research Area

Image Data



Western blot analysis of DGKA in HeLa, A549 lysates using DGKA antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using DGKA antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.