

Product Name: Dnmt1 (10B9) Rabbit Monoclonal Antibody
Catalog #: AMRe10087

Summary

Production Name	Dnmt1 (10B9) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol.
Purification	Affinity purification

Immunogen

Gene Name	DNMT1
Alternative Names	
Gene ID	1786.0
SwissProt ID	P26358.

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	183 kDa


Background

Product Name: Dnmt1 (10B9) Rabbit Monoclonal Antibody
Catalog #: AMRe10087



Research Area

Image Data



Western blot analysis of Dnmt1 in HEK293 , NIH/3T3, C6 lysates using Dnmt1 antibody.

Note

For research use only.