

Product Name: ERp57 (7K11) Rabbit Monoclonal Antibody
Catalog #: AMRe10613

Summary

Production Name	ERp57 (7K11) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	PDIA3
Alternative Names	ERp57; ERp60; ERp61; Glucose Regulated Protein 58 Kd; HsT17083; p58; PDIA3; Phospholipase C alpha; PI PLC;
Gene ID	2923.0
SwissProt ID	P30101.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	57kDa


Product Name: ERp57 (7K11) Rabbit Monoclonal Antibody
Catalog #: AMRe10613

Background

Detected in the flagellum and head region of spermatozoa (at protein level).

Research Area

Image Data



Western blot analysis of ERp57 expression in HepG2 cell lysate.

Note

For research use only.