

**Product Name: Factor X (14U4) Rabbit Monoclonal Antibody**  
**Catalog #: AMRe10787**



## Summary

|                        |  |
|------------------------|--|
| <b>Production Name</b> | Factor X (14U4) Rabbit Monoclonal Antibody |
| <b>Description</b>     | Rabbit Monoclonal Antibody                 |
| <b>Host</b>            | Rabbit                                     |
| <b>Application</b>     | WB,ELISA                                   |
| <b>Reactivity</b>      | Human                                      |

## Performance

|                     |  |
|---------------------|--|
| <b>Conjugation</b>  | Unconjugated   |
| <b>Modification</b> | Unmodified   |
| <b>Isotype</b>      | IgG  |
| <b>Clonality</b>    | Monoclonal   |
| <b>Form</b>         | Liquid   |
| <b>Storage</b>      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.   |
| <b>Buffer</b>       | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle. |
| <b>Purification</b> | Affinity purification  |

## Immunogen

|                          |  |
|--------------------------|--|
| <b>Gene Name</b>         | F10  |
| <b>Alternative Names</b> | Activated factor Xa heavy chain; Coagulation factor X; F10; factor Xa; FXA; Prothrombinase; Stuart factor; |
| <b>Gene ID</b>           | 2159.0   |
| <b>SwissProt ID</b>      | P00742.  |

## Application

|                         |                  |
|-------------------------|------------------|
| <b>Dilution Ratio</b>   | WB 1:1000-1:5000 |
| <b>Molecular Weight</b> | 55kDa            |

**Product Name: Factor X (14U4) Rabbit Monoclonal Antibody**  
**Catalog #: AMRe10787**




## Background

Factor Xa is a vitamin K-dependent glycoprotein that converts prothrombin to thrombin in the presence of factor Va, calcium and phospholipid during blood clotting. Factor Xa is a vitamin K-dependent glycoprotein that converts prothrombin to thrombin in the presence of factor Va, calcium and phospholipid during blood clotting.

## Research Area

## Image Data



Western blot analysis of Factor X expression in HepG2 cell lysate.

## Note

For research use only.