

Product Name: FOXF1 (8B18) Rabbit Monoclonal Antibody
Catalog #: AMRe11086

Summary

Production Name	FOXF1 (8B18) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	FOXF1
Alternative Names	FKHL5; foxf1; FREAC1; HFH 8;
Gene ID	2294.0
SwissProt ID	Q12946.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	40kDa


Product Name: FOXF1 (8B18) Rabbit Monoclonal Antibody
Catalog #: AMRe11086

Background

Probable transcription activator for a number of lung-specific genes. Probable transcription activator for a number of lung-specific genes.

Research Area

Image Data



Western blot analysis of FOXF1 expression in HeLa cell lysate.

Note

For research use only.