

Summary

Production Name	FRA1 (4C9) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	FOSL1
Alternative Names	FOS L1; FOS like antigen 1; FOSL1; FRA 1;
Gene ID	8061.0
SwissProt ID	P15407.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	29kDa

Background

Product Name: FRA1 (4C9) Rabbit Monoclonal Antibody
Catalog #: AMRe11127




Belongs to the bZIP family. Fos subfamily.

Contains 1 bZIP domain.

Research Area

Image Data



Western blot analysis of FRA1 expression in HL60 cell lysate.

Note

For research use only.