

Product Name: FTO (1H6) Rabbit Monoclonal Antibody
Catalog #: AMRe11174



Summary

Production Name	FTO (1H6) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	FTO
Alternative Names	
Gene ID	79068.0
SwissProt ID	Q9C0B1.A synthetic peptide of human FTO

Application

Dilution Ratio	WB: 1:2000-1:10000
Molecular Weight	58kDa


Background

Product Name: FTO (1H6) Rabbit Monoclonal Antibody
Catalog #: AMRe11174



Research Area

Image Data



Western blot detection of FTO in HeLa, Jurkat cell lysates using FTO antibody (1:500 diluted).

Note

For research use only.