

Product Name: GCDFP 15 (8J14) Rabbit Monoclonal Antibody
Catalog #: AMRe11343



Summary

Production Name	GCDFP 15 (8J14) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	PIP
Alternative Names	BRST2; GCDGP15; gp17; GPIP4; PIP; Prolactin induced protein; SABP; Secretory actin binding protein;
Gene ID	5304.0
SwissProt ID	P12273.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	17kDa

Product Name: GCDFP 15 (8J14) Rabbit Monoclonal Antibody
Catalog #: AMRe11343




Background

Belongs to the PIP family.

Research Area

Image Data



Western blot analysis of GCDFP 15 expression in T47D cell lysate.

Note

For research use only.