

**Product Name: GSK3 alpha/beta (3Q8) Rabbit
Monoclonal Antibody
Catalog #: AMRe11819**

Summary

Production Name	GSK3 alpha/beta (3Q8) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	GSK3A
Alternative Names	
Gene ID	2931.0
SwissProt ID	P49840.Recombinant protein of human GSK3 alpha

Application

Dilution Ratio	WB: 1:2000
Molecular Weight	51,47kDa


Background

**Product Name: GSK3 alpha/beta (3Q8) Rabbit
Monoclonal Antibody
Catalog #: AMRe11819**



Research Area

Image Data



Western blot analysis of extracts from A549 cells using RM1200 at 1:1000.

Note

For research use only.