

Product Name: GTPBP4 (14J8) Rabbit Monoclonal Antibody
Catalog #: AMRe11848

Summary

Production Name	GTPBP4 (14J8) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	GTPBP4
Alternative Names	CRFG; GTP binding protein 4; Gtpbp4; NGB; NOG1;
Gene ID	23560.0
SwissProt ID	Q9BZE4.A synthetic peptide of human GTPBP4/NOG1

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	74kDa


Background

Product Name: GTPBP4 (14J8) Rabbit Monoclonal Antibody
Catalog #: AMRe11848

Involved in the biogenesis of the 60S ribosomal subunit. Involved in the biogenesis of the 60S ribosomal subunit.

Research Area

Image Data



Western blot detection of GTPBP4/NOG1 in A549,HL-60 using GTPBP4/NOG1 antibody(1:1000 diluted)

Note

For research use only.