

Product Name: HES5 (17S12) Rabbit Monoclonal Antibody
Catalog #: AMRe11993



Summary

Production Name	HES5 (17S12) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	HES5
Alternative Names	HES 5; Hes5;
Gene ID	388585.0
SwissProt ID	Q5TA89.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	18kDa


Product Name: HES5 (17S12) Rabbit Monoclonal Antibody
Catalog #: AMRe11993

Background

Transcriptional repressor of genes that require a bHLH protein for their transcription. Transcriptional repressor of genes that require a bHLH protein for their transcription. Plays an important role as neurogenesis negative regulator (By similarity).

Research Area

Image Data



Western blot analysis of HES5 expression in (1) 293 cell lysate; (2) Mouse spleen lysate; (3) C6 cell lysate.

Note

For research use only.