

Product Name: Histone H1.2 (19S17) Rabbit Monoclonal Antibody
Catalog #: AMRe12053



Summary

Production Name	Histone H1.2 (19S17) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	H1-2
Alternative Names	
Gene ID	3006.0
SwissProt ID	P16403.A synthetic peptide of human Histone H1.2

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	21kDa


Background

Product Name: Histone H1.2 (19S17) Rabbit Monoclonal Antibody
Catalog #: AMRe12053



Research Area

Image Data



Western blot detection of Histone H1.2 in HeLa,C6,MCF7 cell lysates using Histone H1.2 antibody(1:1000 diluted).

Note

For research use only.