

**Product Name: Histone H2A.X (13P7) Rabbit
Monoclonal Antibody
Catalog #: AMRe12056**

Summary

Production Name	Histone H2A.X (13P7) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	H2AX
Alternative Names	
Gene ID	3014.0
SwissProt ID	P16104.A synthetic peptide of human Histone H2A.X

Application


Dilution Ratio	WB: 1:1000-1:5000
Molecular Weight	15kDa

Background

**Product Name: Histone H2A.X (13P7) Rabbit
Monoclonal Antibody
Catalog #: AMRe12056**

Research Area

Image Data



Western blot detection of Histone H2A.X in 3T3,HeLa cell lysates using Histone H2A.X antibody(1:1000 diluted).

Note

For research use only.