

Product Name: HMGB1 (2P7) Rabbit Monoclonal Antibody
Catalog #: AMRe12109

Summary

Production Name	HMGB1 (2P7) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% New type preservative N and 0.05% BSA
Purification	Affinity purification

Immunogen

Gene Name	HMGB1
Alternative Names	
Gene ID	3146.0
SwissProt ID	P09429.

Application


Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100
Molecular Weight	25 kDa

Product Name: HMGB1 (2P7) Rabbit Monoclonal Antibody
Catalog #: AMRe12109


Background

Research Area

Image Data



Western blot analysis of HMGB1 in K562, rat Brain, C6, 3T3, HeLa lysates using HMGB1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using HMGB1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.