

Product Name: hnRNP G (8L8) Rabbit Monoclonal Antibody
Catalog #: AMRe12144



Summary

Production Name	hnRNP G (8L8) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	RBMX
Alternative Names	
Gene ID	27316.0
SwissProt ID	P38159.A synthetic peptide of human hnRNP G

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	42kDa


Background

Product Name: hnRNP G (8L8) Rabbit Monoclonal Antibody
Catalog #: AMRe12144



Research Area

Image Data



Western blot detection of hnRNP G in K562, HeLa cell lysates using hnRNP G antibody (1:1000 diluted).

Note

For research use only.