

Product Name: Huntingtin (2P4) Rabbit Monoclonal Antibody
Catalog #: AMRe12284



Summary

Production Name	Huntingtin (2P4) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	HTT
Alternative Names	HD; HDH; HTT; Huntingtin; IT15; ZHD;
Gene ID	3064.0
SwissProt ID	P42858.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	348kDa

Product Name: Huntingtin (2P4) Rabbit Monoclonal Antibody
Catalog #: AMRe12284




Background

May play a role in microtubule-mediated transport or vesicle function. [Huntingtin]: May play a role in microtubule-mediated transport or vesicle function.

Research Area

Image Data



Western blot analysis of Huntingtin expression in HeLa cell lysate.

Note

For research use only.