

Product Name: IGF1 Receptor (18P7) Rabbit Monoclonal Antibody
Catalog #: AMRe12421



Summary

Production Name	IGF1 Receptor (18P7) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	IGF1R
Alternative Names	
Gene ID	3480.0
SwissProt ID	P08069.Recombinant protein of human IGF1 Receptor

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	155kDa


Background

Product Name: IGF1 Receptor (18P7) Rabbit Monoclonal Antibody
Catalog #: AMRe12421



Research Area

Image Data



Western blot detection of IGF1 Receptor in HeLa,A549,MCF-7,U251 cell lysates using IGF1 Receptor antibody(1:1000 diluted).

Note

For research use only.