

**Product Name: IL18 binding protein (18K17) Rabbit
Monoclonal Antibody
Catalog #: AMRe12513**



Summary

Production Name	IL18 binding protein (18K17) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	IL18BP
Alternative Names	IL18 binding protein; IL18 Bpa; IL18BP; IL18BPa; Interleukin 18 binding protein; MC51L 53L 54L homolog gene product; Tadekinig alfa;
Gene ID	10068.0
SwissProt ID	O95998.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	21kDa

**Product Name: IL18 binding protein (18K17) Rabbit
Monoclonal Antibody
Catalog #: AMRe12513**




Background

Isoform A binds to IL-18 and inhibits its activity. Functions as an inhibitor of the early TH1 cytokine response. Isoform A binds to IL-18 and inhibits its activity. Functions as an inhibitor of the early TH1 cytokine response.

Research Area

Image Data



Western blot analysis of IL18 binding protein in HeLa cell lysate.

Note

For research use only.