

Summary

Production Name	IL9 (14B18) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	IL9
Alternative Names	HP40; IL9; Interleukin 9 MCGF; P40; TCGF 3;
Gene ID	3578.0
SwissProt ID	P15248.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	16kDa

Background


Product Name: IL9 (14B18) Rabbit Monoclonal Antibody
Catalog #: AMRe12574



Supports IL-2 independent and IL-4 independent growth of helper T-cells. Supports IL-2 independent and IL-4 independent growth of helper T-cells.

Research Area

Image Data



Western blot analysis of IL9 in Recombinant Human IL9 protein cell lysate.

Note

For research use only.