

**Product Name: Intra Acrosomal Protein (7H17) Rabbit
Monoclonal Antibody
Catalog #: AMRe12688**



Summary

Production Name	Intra Acrosomal Protein (7H17) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	ACRV1
Alternative Names	ACRV1; Msa63; Sp10; SPACA2;
Gene ID	56.0
SwissProt ID	P26436.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	28kDa

**Product Name: Intra Acrosomal Protein (7H17) Rabbit
Monoclonal Antibody
Catalog #: AMRe12688**




Background

N/A

Research Area

Image Data



Western blot analysis of Intra Acrosomal Protein expression in Human testis lysate.

Note

For research use only.