

Product Name: IQGAP2 (18D18) Rabbit Monoclonal Antibody
Catalog #: AMRe12721



Summary

Production Name	IQGAP2 (18D18) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	IQGAP2
Alternative Names	lqqap2;
Gene ID	10788.0
SwissProt ID	Q13576.Recombinant protein of human IQGAP2

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	181kDa

Background


Product Name: IQGAP2 (18D18) Rabbit Monoclonal Antibody
Catalog #: AMRe12721



Binds to activated CDC42 and RAC1 but does not seem to stimulate their GTPase activity. Associates with calmodulin. Binds to activated CDC42 and RAC1 but does not seem to stimulate their GTPase activity. Associates with calmodulin.

Research Area

Image Data



Western blot detection of IQGAP2 in A549,HL-60 using IQGAP2 antibody(1:1000 diluted)

Note

For research use only.