

Product Name: IQGAP2 (18P19) Rabbit Monoclonal Antibody
Catalog #: AMRe12722



Summary

Production Name	IQGAP2 (18P19) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	IQGAP2
Alternative Names	lqqap2;
Gene ID	10788.0
SwissProt ID	Q13576.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	181kDa

Product Name: IQGAP2 (18P19) Rabbit Monoclonal Antibody
Catalog #: AMRe12722




Background

Binds to activated CDC42 and RAC1 but does not seem to stimulate their GTPase activity. Associates with calmodulin. Binds to activated CDC42 and RAC1 but does not seem to stimulate their GTPase activity. Associates with calmodulin.

Research Area

Image Data



Western blot analysis of IQGAP2 expression in HepG2 cell lysate.

Note

For research use only.