

Product Name: IRAK2 (2Q15) Rabbit Monoclonal Antibody
Catalog #: AMRe12726

Summary

Production Name	IRAK2 (2Q15) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	IRAK2
Alternative Names	Interleukin 1 receptor associated kinase 2; IRAK 2; Irak2;
Gene ID	3656.0
SwissProt ID	O43187.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	69kDa


Product Name: IRAK2 (2Q15) Rabbit Monoclonal Antibody
Catalog #: AMRe12726

Background

Binds to the IL-1 type I receptor following IL-1 engagement, triggering intracellular signaling cascades leading to transcriptional up-regulation and mRNA stabilization. Binds to the IL-1 type I receptor following IL-1 engagement, triggering intracellular signaling cascades leading to transcriptional up-regulation and mRNA stabilization.

Research Area

Image Data



Western blot analysis of IRAK2 expression in HeLa cell lysate.

Note

For research use only.