

Product Name: LDHA (4H19) Rabbit Monoclonal Antibody
Catalog #: AMRe13262

Summary

Production Name	LDHA (4H19) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	LDHA
Alternative Names	LDHA; LDH-M; LDH1; PIG19;
Gene ID	3939.0
SwissProt ID	P00338.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	37kDa


Product Name: LDHA (4H19) Rabbit Monoclonal Antibody
Catalog #: AMRe13262

Background

Fermentation; pyruvate fermentation to lactate; (S)-lactate from pyruvate: step 1/1.

Research Area

Image Data



Western blot analysis of LDHA expression in MCF7 cell lysate.

Note

For research use only.