

Product Name: LPCAT1 (13M5) Rabbit Monoclonal Antibody
Catalog #: AMRe13387



Summary

Production Name	LPCAT1 (13M5) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	LPCAT1
Alternative Names	
Gene ID	79888.0
SwissProt ID	Q8NF37.Recombinant protein of human LPCAT1

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	59kDa


Background

Product Name: LPCAT1 (13M5) Rabbit Monoclonal Antibody
Catalog #: AMRe13387



Research Area

Image Data



Western blot detection of LPCAT1 in A549 cell lysates using LPCAT1 antibody(1:1000 diluted).

Note

For research use only.