

Product Name: LRG1 (8F15) Rabbit Monoclonal Antibody
Catalog #: AMRe13418

Summary

Production Name	LRG1 (8F15) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	LRG1
Alternative Names	A2GL; HMFT1766; LRG; LRG1;
Gene ID	116844.0
SwissProt ID	P02750.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	38kDa


Product Name: LRG1 (8F15) Rabbit Monoclonal Antibody
Catalog #: AMRe13418

Background

Transforming growth factor beta receptor binding.

Research Area

Image Data



Western blot analysis of LRG1 expression in Human plasma lysate.

Note

For research use only.