

Product Name: MAS1L (13G3) Rabbit Monoclonal Antibody
Catalog #: AMRe13652

Summary

Production Name	MAS1L (13G3) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	MAS1L
Alternative Names	MAS L; MAS R; MAS1L; MRG;
Gene ID	116511.0
SwissProt ID	P35410.A synthetic peptide of human MAS1L

Application

Dilution Ratio	WB: 1:2000-1:10000
Molecular Weight	42kDa


Background

Product Name: MAS1L (13G3) Rabbit Monoclonal Antibody
Catalog #: AMRe13652

Belongs to the G-protein coupled receptor 1 family. Mas subfamily.

Research Area

Image Data



Western blot detection of MAS1L in K562,Rat Brain,C6,3T3,HeLa cell lysates using MAS1L antibody(1:1000 diluted).

Note

For research use only.