

Product Name: MLK3 (1E14) Rabbit Monoclonal Antibody
Catalog #: AMRe13954

Summary

Production Name	MLK3 (1E14) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% New type preservative N and 0.05% BSA
Purification	Affinity purification

Immunogen

Gene Name	MAP3K11
Alternative Names	
Gene ID	4296.0
SwissProt ID	Q16584.

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	93 kDa


Background

Product Name: MLK3 (1E14) Rabbit Monoclonal Antibody
Catalog #: AMRe13954



Research Area

Image Data



Western blot analysis of MLK3 in HepG2, A431 lysates using MLK3 antibody.

Note

For research use only.