

Product Name: MSRB3 (4D16) Rabbit Monoclonal Antibody
Catalog #: AMRe14184

Summary

Production Name	MSRB3 (4D16) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	MSRB3
Alternative Names	DFNB74; MsrB3;
Gene ID	253827.0
SwissProt ID	Q8IXL7.Recombinant protein of human MSRB3

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	21kDa


Background

Product Name: MSRB3 (4D16) Rabbit Monoclonal Antibody
Catalog #: AMRe14184

Catalyzes the reduction of free and protein-bound methionine sulfoxide to methionine. Isoform 2 is essential for hearing.
Catalyzes the reduction of free and protein-bound methionine sulfoxide to methionine. Isoform 2 is essential for hearing.

Research Area

Image Data



Western blot detection of MSRB3 in C6 cell lysates using MSRB3 antibody(1:1000 diluted).

Note

For research use only.