

Product Name: MTX3 (10I6) Rabbit Monoclonal Antibody
Catalog #: AMRe14224

Summary

Production Name	MTX3 (10I6) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	MTX3
Alternative Names	
Gene ID	345778.0
SwissProt ID	Q5HYI7.A synthetic peptide of human MTX3

Application

Dilution Ratio	WB: 1:1000-1:5000
Molecular Weight	35kDa


Background

Product Name: MTX3 (10I6) Rabbit Monoclonal Antibody
Catalog #: AMRe14224



Research Area

Image Data



Western blot detection of MTX3 in mouse brain cell lysates using MTX3 antibody(1:1000 diluted).

Note

For research use only.