

**Product Name: Myeloperoxidase (7118) Rabbit
Monoclonal Antibody
Catalog #: AMRe14281**

Summary

Production Name	Myeloperoxidase (7118) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% New type preservative N and 0.05% BSA
Purification	Affinity purification

Immunogen

Gene Name	MPO
Alternative Names	
Gene ID	4353.0
SwissProt ID	P05164.

Application


Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100
Molecular Weight	84 kDa

**Product Name: Myeloperoxidase (7I18) Rabbit
Monoclonal Antibody
Catalog #: AMRe14281**


Background

Research Area

Image Data



Western blot analysis of Myeloperoxidase in 3T3 lysates using Myeloperoxidase antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Myeloperoxidase antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.