

Product Name: MYH11 (10S12) Rabbit Monoclonal Antibody
Catalog #: AMRe14289



Summary

Production Name	MYH11 (10S12) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	MYH11
Alternative Names	AAT4; FAA4; MYH11; SM1; SM2; SMHC; SMMHC;
Gene ID	4629.0
SwissProt ID	P35749.

Application

Dilution Ratio	WB 1:1000-1:5000
Molecular Weight	227kDa

Product Name: MYH11 (10S12) Rabbit Monoclonal Antibody
Catalog #: AMRe14289




Background

Muscle contraction. Muscle contraction.

Research Area

Image Data



Western blot analysis of Myosin, smooth muscle heavy chain 1 and 2 expression in Human bladder lysate.

Note

For research use only.