

Product Name: MYH11 (6C13) Rabbit Monoclonal Antibody
Catalog #: AMRe14290



Summary

Production Name	MYH11 (6C13) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	MYH11
Alternative Names	AAT4; FAA4; MYH11; Myosin 11; Myosin heavy chain 11; Myosin heavy chain smooth muscle isoform; Myosin heavy polypeptide 11 smooth muscle; SMHC; SMMHC;
Gene ID	4629.0
SwissProt ID	P35749.A synthetic peptide of human Myosin, smooth muscle heavy chain 1 and 2

Application

Dilution Ratio	WB: 1:2000-1:10000
Molecular Weight	227kDa


Product Name: MYH11 (6C13) Rabbit Monoclonal Antibody
Catalog #: AMRe14290

Background

Muscle contraction. Muscle contraction.

Research Area

Image Data



Western blot analysis of extracts from Mouse spleen tissue using RM5298 at 1:1000.

Note

For research use only.