

Product Name: NeuN (11E16) Rabbit Monoclonal Antibody
Catalog #: AMRe14584

Summary

Production Name	NeuN (11E16) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	RBFOX3
Alternative Names	
Gene ID	146713.0
SwissProt ID	A6NFN3.A synthetic peptide of human NeuN

Application


Dilution Ratio	WB: 1:1000
Molecular Weight	34kDa

Background

Product Name: NeuN (11E16) Rabbit Monoclonal Antibody
Catalog #: AMRe14584

Research Area

Image Data



Western blot detection of NeuN in Rat Brain lysates using NeuN antibody(1:1000 diluted).

Note

For research use only.