

Product Name: Oct-2 (16S17) Rabbit Monoclonal Antibody
Catalog #: AMRe15101

Summary

Production Name	Oct-2 (16S17) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	POU2F2
Alternative Names	
Gene ID	5452.0
SwissProt ID	P09086.Recombinant protein of human Oct-2

Application


Dilution Ratio	WB: 1:1000
Molecular Weight	51kDa

Background

Product Name: Oct-2 (16S17) Rabbit Monoclonal Antibody
Catalog #: AMRe15101

Research Area

Image Data



Western blot analysis of extracts from Jurkat cells using RM1245 at 1:1000.

Note

For research use only.