

Product Name: OTUB1 (7P2) Rabbit Monoclonal Antibody
Catalog #: AMRe15534



Summary

Production Name	OTUB1 (7P2) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	OTUB1
Alternative Names	
Gene ID	55611.0
SwissProt ID	Q96FW1.A synthetic peptide of human OTUB1

Application


Dilution Ratio	WB: 1:2000
Molecular Weight	31kDa

Background

Product Name: OTUB1 (7P2) Rabbit Monoclonal Antibody
Catalog #: AMRe15534

Research Area

Image Data



Western blot detection of OTUB1 in K562,Rat Brain,C6,3T3,HeLa cell lysates using OTUB1 antibody(1:1000 diluted).

Note

For research use only.