

Product Name: p27 Kip 1 (18B10) Rabbit Monoclonal Antibody
Catalog #: AMRe15592



Summary

Production Name	p27 Kip 1 (18B10) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	CDKN1B
Alternative Names	
Gene ID	1027.0
SwissProt ID	P46527.A synthetic peptide of human p27 KIP 1

Application

Dilution Ratio	WB: 1:1000-1:10000
Molecular Weight	22kDa


Background

Product Name: p27 Kip 1 (18B10) Rabbit Monoclonal Antibody
Catalog #: AMRe15592



Research Area

Image Data



Western blot detection of p27 KIP 1 in 3T3,HeLa cell lysates using p27 KIP 1 antibody(1:1000 diluted).

Note

For research use only.