

Product Name: P4HB (15Y14) Rabbit Monoclonal Antibody
Catalog #: AMRe15635



Summary

Production Name	P4HB (15Y14) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	P4HB
Alternative Names	
Gene ID	5034.0
SwissProt ID	P07237.A synthetic peptide of human P4HB

Application

Dilution Ratio	WB: 1:1000-1:5000
Molecular Weight	57kDa


Background

Product Name: P4HB (15Y14) Rabbit Monoclonal Antibody
Catalog #: AMRe15635



Research Area

Image Data



Western blot detection of P4HB in HCT116, HeLa, HepG2, Jurkat, MCF-7 cell lysates using P4HB antibody (1:1000 diluted).

Note

For research use only.