

Product Name: p62/SQSTM1 (12N5) Rabbit Monoclonal Antibody
Catalog #: AMRe15658

Summary

Production Name	p62/SQSTM1 (12N5) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	SQSTM1
Alternative Names	OSIL; Oxidative stress induced like; p60; p62; p62B; Paget disease of bone 3;PDB 3; PDB3;
Gene ID	8878.0
SwissProt ID	Q13501.

Application

Dilution Ratio	WB 1:5000~1:10000
Molecular Weight	48kDa


Background

Autophagy receptor that interacts directly with both the cargo to become degraded and an autophagy modifier of the MAP1 LC3 family. Required both for the formation and autophagic degradation of polyubiquitin-containing bodies, called ALIS (aggresome-like induced structures) and links ALIS to the autophagic machinery. Involved in midbody ring degradation. May regulate the activation of NFKB1 by TNF-alpha, nerve growth factor (NGF) and interleukin-1. Autophagy receptor required for selective macroautophagy (aggrephagy). Functions as a bridge between polyubiquitinated cargo and autophagosomes. Interacts directly with both the cargo to become degraded and an autophagy modifier of the MAP1 LC3 family (PubMed: [16286508](http://www.uniprot.org/citations/16286508), PubMed: [20168092](http://www.uniprot.org/citations/20168092), PubMed: [24128730](http://www.uniprot.org/citations/24128730), PubMed: [28404643](http://www.uniprot.org/citations/28404643), PubMed: [22622177](http://www.uniprot.org/citations/22622177)). Along with WDFY3, involved in the formation and autophagic degradation of cytoplasmic ubiquitin-containing inclusions (p62 bodies, ALIS/aggresome-like induced structures). Along with WDFY3, required to recruit ubiquitinated proteins to PML bodies in the nucleus (PubMed: [24128730](http://www.uniprot.org/citations/24128730), PubMed: [20168092](http://www.uniprot.org/citations/20168092)). May regulate the activation of NFKB1 by TNF-alpha, nerve growth factor (NGF) and interleukin-1. May play a role in titin/TTN downstream signaling in muscle cells. May regulate signaling cascades through ubiquitination. Adapter that mediates the interaction between TRAF6 and CYLD (By similarity). May be involved in cell differentiation, apoptosis, immune response and regulation of K(+) channels. Involved in endosome organization by retaining vesicles in the perinuclear cloud: following ubiquitination by RNF26, attracts specific vesicle-associated adapters, forming a molecular bridge that restrains cognate vesicles in the perinuclear region and organizes the endosomal pathway for efficient cargo transport (PubMed: [27368102](http://www.uniprot.org/citations/27368102)). Promotes relocalization of 'Lys-63'-linked ubiquitinated STING1 to autophagosomes (PubMed: [29496741](http://www.uniprot.org/citations/29496741)). Acts as an activator of the NFE2L2/NRF2 pathway via interaction with KEAP1: interaction inactivates the BCR(KEAP1) complex, promoting nuclear accumulation of NFE2L2/NRF2 and subsequent expression of cytoprotective genes (PubMed: [20452972](http://www.uniprot.org/citations/20452972), PubMed: [28380357](http://www.uniprot.org/citations/28380357)).

Research Area

Image Data

Product Name: p62/SQSTM1 (12N5) Rabbit Monoclonal Antibody
Catalog #: AMRe15658



Western blot analysis of p62/SQSTM1 expression in SKBR-3 cell lysate.

Note

For research use only.