

Product Name: PD L2 (12P7) Rabbit Monoclonal Antibody
Catalog #: AMRe15873

Summary

Production Name	PD L2 (12P7) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	PDCD1LG2
Alternative Names	
Gene ID	80380.0
SwissProt ID	Q9BQ51.A synthetic peptide of human PD-L2

Application

Dilution Ratio	WB: 1:1000-1:5000
Molecular Weight	31kDa


Background

Product Name: PD L2 (12P7) Rabbit Monoclonal Antibody
Catalog #: AMRe15873



Research Area

Image Data



Western blot detection of PD-L2 in K562,C6 cell lysates using PD-L2 antibody(1:1000 diluted).

Note

For research use only.