

Product Name: Pentraxin 3 (14E14) Rabbit Monoclonal Antibody
Catalog #: AMRe15962



Summary

Production Name	Pentraxin 3 (14E14) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	PTX3
Alternative Names	Pentaxin 3; PTX3; TNFAIP5; TSG14;
Gene ID	5806.0
SwissProt ID	P26022.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	42kDa

Product Name: Pentraxin 3 (14E14) Rabbit Monoclonal Antibody
Catalog #: AMRe15962




Background

Plays a role in the regulation of innate resistance to pathogens, inflammatory reactions, possibly clearance of self-components and female fertility. Plays a role in the regulation of innate resistance to pathogens, inflammatory reactions, possibly clearance of self-components and female fertility.

Research Area

Image Data



Western blot analysis of Pentraxin 3 expression in HUVEC cell lysate.

Note

For research use only.