

Product Name: PGP9.5 (5O2) Rabbit Monoclonal Antibody
Catalog #: AMRe16043

Summary

Production Name	PGP9.5 (5O2) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	UCHL1
Alternative Names	UCH-L1; PGP9.5; UCHL1;
Gene ID	7345.0
SwissProt ID	P09936.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	25kDa


Product Name: PGP9.5 (5O2) Rabbit Monoclonal Antibody
Catalog #: AMRe16043

Background

Involved both in the processing of ubiquitin precursors and of ubiquitinated proteins. This enzyme is a thiol protease that recognizes and hydrolyzes a peptide bond at the C-terminal glycine of ubiquitin. Ubiquitin-protein hydrolase involved both in the processing of ubiquitin precursors and of ubiquitinated proteins (Probable). This enzyme is a thiol protease that recognizes and hydrolyzes a peptide bond at the C-terminal glycine of ubiquitin (PubMed:9774100, PubMed:8639624, PubMed:12408865, PubMed:23359680). Also binds to free monoubiquitin and may prevent its degradation in lysosomes (By similarity). The homodimer may have ATP-independent ubiquitin ligase activity (PubMed:12408865).

Research Area

Image Data



Western blot analysis of PGP9.5 expression in SH-SY5Y cell lysate.

Note

For research use only.