

Product Name: Pleckstrin (2W13) Rabbit Monoclonal Antibody
Catalog #: AMRe16262



Summary

Production Name	Pleckstrin (2W13) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	PLEK
Alternative Names	p47; Pleckstrin; PLEK;
Gene ID	5341.0
SwissProt ID	P08567.A synthetic peptide of human Pleckstrin

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	40kDa

Background


Product Name: Pleckstrin (2W13) Rabbit Monoclonal Antibody
Catalog #: AMRe16262



Major protein kinase C substrate of platelets. Major protein kinase C substrate of platelets.

Research Area

Image Data



Western blot detection of Pleckstrin in HeLa,HL-60 cell lysates using Pleckstrin antibody(1:1000 diluted).

Note

For research use only.