

Product Name: PMP22 (2M9) Rabbit Monoclonal Antibody
Catalog #: AMRe16308

Summary

Production Name	PMP22 (2M9) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	PMP22
Alternative Names	CMT1A; CMT1E; DSS; GAS3; Growth Arrest Specific 3; HMSNIA; HNPP; Peripheral myelin protein 22; PMP22; Sp110; Trembler;
Gene ID	5376.0
SwissProt ID	Q01453.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	18kDa


Product Name: PMP22 (2M9) Rabbit Monoclonal Antibody
Catalog #: AMRe16308

Background

Involved in growth regulation, and in myelination in the peripheral nervous system. Might be involved in growth regulation, and in myelination in the peripheral nervous system.

Research Area

Image Data



Western blot analysis of PMP22 expression in 293T cell lysate.

Note

For research use only.