

**Product Name: PPAR gamma (19M16) Rabbit
Monoclonal Antibody
Catalog #: AMRe16409**

Summary

Production Name	PPAR gamma (19M16) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	PPARG
Alternative Names	CIMT1; GLM1; NR1C3; PPAR gamma; PARG; PPARG1; PPARG2;
Gene ID	5468.0
SwissProt ID	P37231.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	58kDa


**Product Name: PPAR gamma (19M16) Rabbit
Monoclonal Antibody
Catalog #: AMRe16409**

Background

Receptor that binds peroxisome proliferators such as hypolipidemic drugs and fatty acids. Once activated by a ligand, the receptor binds to a promoter element in the gene for acyl-CoA oxidase and activates its transcription. It therefore controls the peroxisomal beta-oxidation pathway of fatty acids. Key regulator of adipocyte differentiation and glucose homeostasis. Nuclear receptor that binds peroxisome proliferators such as hypolipidemic drugs and fatty acids. Once activated by a ligand, the nuclear receptor binds to DNA specific PPAR response elements (PPRE) and modulates the transcription of its target genes, such as acyl-CoA oxidase. It therefore controls the peroxisomal beta-oxidation pathway of fatty acids. Key regulator of adipocyte differentiation and glucose homeostasis. ARF6 acts as a key regulator of the tissue-specific adipocyte P2 (aP2) enhancer. Acts as a critical regulator of gut homeostasis by suppressing NF-kappa-B-mediated proinflammatory responses. Plays a role in the regulation of cardiovascular circadian rhythms by regulating the transcription of ARNTL/BMAL1 in the blood vessels (By similarity).

Research Area

Image Data



Western blot analysis of PPAR gamma expression in (1) HeLa cell lysate; (2) PC-3 cell lysate; (3) THP-1 cell lysate.

Note

For research use only.