

Product Name: Prostein (16Q10) Rabbit Monoclonal Antibody
Catalog #: AMRe16529



Summary

Production Name	Prostein (16Q10) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	SLC45A3
Alternative Names	IPCA6; PCANAP2; PCANAP6; PCANAP8; Prostein; PRST; Slc45a3;
Gene ID	85414.0
SwissProt ID	Q96JT2.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	59kDa

Product Name: Prostein (16Q10) Rabbit Monoclonal Antibody
Catalog #: AMRe16529




Background

Marker for prostate cells. May be used, in case of prostate cancers, as a target antigen for prostate carcinomas-directed cytotoxic T-cell lymphocytes.

Research Area

Image Data



Western blot analysis of Prostein expression in LnCaP cell lysate.

Note

For research use only.